

Detection and quantification of water-based aerosols using active open-path FTIR

O. Kira*, R. Linker and Y. Dubowski

Faculty of Civil and Environmental Engineering, Technion-Israel Institute of Technology, Technion City, Haifa 3200003, Israel.

Keywords: Active remote sensing, Thermal infrared, hydrosols.

Presenting author email: ozkira@campus.technion.ac.il

Exposure to airborne particles has an adverse influence on human health. One of the great challenges in aerosol real-time measurements, besides quantification, is to identify their composition. This study suggests using active Open Path Fourier Transfer Infra-Red (OP-FTIR) spectroscopy for quantifying water load in spray cloud as well as solutes load in it.

OP-FTIR operates in the mid-IR range where many organic and inorganic compounds present unique spectral “fingerprints” and hence has great potential for obtaining chemical information of the aerosols. Although a few studies have shown the potential of OP-FTIR for aerosols detection, e.g. measurements of water droplets (Hashmonay & Yost 1999; Wu et al. 2007), It has not been used for investigation of aerosols chemical composition.

Experimental: To minimize wind interference the aerosol cloud was generated within a 20m long polyethylene tunnel open at both ends. A custom agricultural system was used for the dispersion of water-based droplets. In addition to the spectral measurements, off-line measurements of water load in the line of site (LOS) and concentration of the dissolved components in the airborne droplets were conducted. (see figure 1)




Figure 1: Illustration of the experimental setup inside the spray tunnel.

Three sets of spraying experiments generated different hydrosols clouds: (1) tap water only (with varying loads using different nuzzles), (2) aqueous ammonium sulfate solution (0.35%-3.6% wt) and (3) aqueous ethylene glycol solution (0.47%-2.38% wt).

Results

Experiment 1- Water clouds: The wavelength-dependent extinction coefficient of the droplets was found to be similar for all tested nozzles, resulting a linear relationship between spectrum baseline and water load in the LOS.

Experiment 2 - Cloud of aqueous ammonium sulfate droplets: The ammonium sulfate spectral signature was clearly visible in the signal and a quantitative model was obtained (figure 2).




Figure 2: Water-subtracted signals of the ammonium sulfate clouds.

Experiment 3: Cloud of aqueous ethylene glycol droplets. For the semi-volatile ethylene glycol, present in both gas and condense phases, quantification was much more complex and two approaches were developed: (1) using the linear relationship from the first experiment (determination error of 8%), and (2) inverse modeling (determination error of 57%).

Conclusion

This work demonstrates the potential of the OP-FTIR for detecting clouds of water-based aerosols and for quantifying water droplets and solutes present in them at relatively low concentrations.

This work was done under and supported by the Technion Center of Excellence in Exposure Science and Environmental Health (TCEEH). Additionally, this work was supported by the Israel Science Foundation (809/12) and the Environmental Health Foundation (864942). We would like to thank the American Technion Society for their financial support through the generous donation of the Hittman Family.

Hashmonay, R.A. and Yost, M.G. (1999). Environmental Science & Technology, 33(7), 1141–1144.

Wu, C.-F. et al., (2007). Journal of Environmental Science and Health, Part A, 42(8), 1131–1140.